


Posey Reusable Double-Padded Mitts (Mittens) – One Size RESTRAINT

	<p>Indications:</p> <ul style="list-style-type: none"> • Used as a hand restraint for patients who are picking, pulling, scratching or peeling where it exacerbates skin condition, causes self-injury, or disrupts wound site integrity or disrupts life sustaining medical treatment • Depending on degree of restlessness and agitation, may be a good first option
<p>Supply Details</p> <p><u>Mittens</u></p> <p>ITEM ID: 00011537</p> <p>VENDOR ITEM ID: Bowers/Posey 2819</p> <p>Manufacturer web site (video and other information)</p> <p>http://www.tidiproducs.com/product-listing/posey-mitts</p>	<p>Considerations:</p> <ul style="list-style-type: none"> • Mitten is a less restrictive device that should be considered before more restrictive devices are used • Mitten is universal size, padded on both sides, has no finger separators and has a Velcro closure • Machine washable between patient use • Do NOT use with any kind of restraining strap, if the patient requires control of the limb and not just the hand, then an approved limb restraint device must be used • Mitten must be removed according to the standard to assess skin condition and circulation • DO NOT <ul style="list-style-type: none"> ○ Place over watches or jewellery ○ Use if an IV site could be compromised ○ Use on patients with dislocations, fractures or open wounds on the affected limb ○ Apply the mitt too tightly or too loosely, you should be able to insert one finger under the mitten.