

Posey "Twice as Tough" Cuff with D-Ring Non-Locking RESTRAINT Emergency Departments ONLY

	Indications: <ul style="list-style-type: none"> For use in Emergency Department for combative, agitated patients; to prevent disruption of life saving treatment or injury. BLUE cuff for wrist RED cuff for ankle (larger size)
	Considerations: <ul style="list-style-type: none"> Remove watches or jewellery before applying limb restraints Side lying is the recommended positioning for a patient in limb restraints Prone positioning puts the patient at risk for aspiration Do NOT use limb restraints on patients with dislocations, fractures, open wounds, or an IV site on the affected limb Ensure the restraint is not applied too tightly or too loosely. You should be able to insert one finger under the ankle or wrist restraint Attach restraints to bed frame only using the D-ring.
	
Supply Details Blue wrist cuff with D-Rings ITEM ID: 00032154 VENDOR ITEM ID: Bowers/Posey 2790 Red ankle cuff with D-Rings ITEM ID: 00001326 VENDOR ITEM ID: Bowers/Posey 2791	Manufacturer web site (instructions and video)