



## RESPIRATORY SERVICES

DATE CREATED: April 2006

DATE REVIEWED/REVISED: September 2015

## PROCEDURE

TITLE: Pulmonary Diagnostics: Forced Vital Capacity (FVC) (Respiratory Therapy)

NUMBER: B-00-12-12112

RELATED DOCUMENTS:

This material has been prepared solely for use at Providence Health Care (PHC), Provincial Health Services Authority (PHSA) and Vancouver Coastal Health (VCH). PHC, PHSA and VCH accept no responsibility for use of this material by any person or organization not associated with PHC, PHSA and VCH. A printed copy of this document may not reflect the current electronic version.

### SITE APPLICABILITY:

ST. PAUL'S HOSPITAL  
MOUNT SAINT JOSEPH HOSPITAL

### GENERAL INFORMATION:

Forced Vital Capacity (FVC) is the maximal volume of air exhaled with maximally forced effort from a maximal inspiration. FVC measurement is the first measurement in the *Complete Pulmonary Function* testing sequence.

### INDICATIONS:

- To determine the presence or absence of lung disease
- To measure the severity of disease on lung function
- To quantify the effects of exposure to occupational or environmental triggers
- To see if therapy is useful or if there are any side effects
- To evaluate the risk for surgical procedures
- To measure the level of disability or impairment

### CONTRAINDICATIONS:

- Hemoptysis of unknown origin
- Pneumothorax
- Unstable cardiovascular status (recent MI or PE)
- Thoracic, abdominal, or cerebral aneurysms
- Recent thoracic, abdominal, or ocular surgeries
- Acute nausea or vomiting

### CAUTIONS:

Patients should be observed for signs of dizziness and trials stopped accordingly. Instruct patient to not tense their shoulders and neck in order to reduce the likelihood of syncope.

## SPECIAL CONSIDERATIONS:

If patients require supplemental O<sub>2</sub> they should have access to wall O<sub>2</sub> in order to conserve their own supply. Patients with significant kyphosis or scoliosis may have difficulty with mouthpiece height and angle.

## REQUIRED SUPPLIES & EQUIPMENT:

- Jaeger Masterscreen Spirometer
- Microgard Filter
- Silicone Mouthpiece
- Nose clip

## PROCEDURE:

1. Check to make sure the patient does not have any contraindications to testing.
2. The patient should be sitting comfortably, with both feet flat on the floor and the mouthpiece adjusted to a comfortable height.
3. Explain procedure to patient and demonstrate appropriate technique.
4. Double click on the **Spirometry** icon on the main screen.
5. Press **F3** to start measurement.
6. Check to make sure that there is no drift in the baseline prior to placing patient on the mouthpiece.
7. Ask the patient to place the mouthpiece in their mouth and place the nose clips on their nose.
8. Instruct the patient to breathe normally in and out through their mouth.
9. Have the patient take a large breath in and without hesitation exhale forcefully and completely (minimum 6 seconds, maximum 15 seconds) then quickly breathe in until they are full.
10. Press **F7** when maneuver is completed to calculate results.
11. Give feedback and encouragement when appropriate.
12. Watch patient for signs of distress during maneuver – stop if patient becomes dizzy as syncope could follow.
13. A minimum of 3 acceptable maneuvers is required (maximum of 8).
14. Repeatability is achieved when the difference between the two best FVC and FEV<sub>1</sub> is less than 0.15 L and the peak flow is within 10%.
15. The largest combined value of FVC and FEV<sub>1</sub> from acceptable maneuvers should be reported
16. Press **F12** to save results and exit screen.

## REFERENCES:

1. American Thoracic Society Series *ATS/ERS Task Force: Standardization of Spirometry*. European Respiratory Journal 2005 Volume 26
2. *ATS Pulmonary Function Laboratory Management and Procedure Manual (2005)*.

## REVIEWED BY:

1. Respiratory Therapist, Pulmonary Diagnostics, PHC
2. Pulmonary Diagnostics Coordinator, Respiratory Services, PHC
3. Medical Director, Pulmonary Diagnostics, PHC