

Wrist and Ankle Cuffs: Reusable Non- Locking with Quick Release Buckles Restraint

<p>Wrist</p> 	<p>Indications:</p> <p>For use with behaviours where it is necessary to obtain or maintain control of one or more limbs of patients who continually disrupt lifesaving treatment (e.g., chronic tube pulling which breaks the integrity of an IV line or self-extubation)</p>
<p>Ankle</p> 	<p>Considerations:</p> <ul style="list-style-type: none"> • Remove watches or jewellery before applying limb restraints • Side-Lying is the preferred positioning when a patient is in limb restraints. Prone positioning places patient at risk for aspiration • Attachment to the bed frame must be made using a quick release buckle as recommended by the manufacturer • DO NOT use limb restraints on patients with dislocations, fractures, open wounds or an IV site on the affected limb • Ensure the restraint is not applied too tightly or too loosely. You should be able to insert one finger under the ankle or wrist restraint • Attach restraints to bed frame only
<p>Supply Details</p> <p>Blue wrist cuff ITEM ID: 00090260 VENDOR ITEM ID: Bowers/Posey 2790Q</p> <p>Red ankle cuff ITEM ID: 00090261 VENDOR ITEM ID: Bowers/Posey 2791Q</p>	<p>Manufacturers web site</p>

This material has been prepared solely for use at Providence Health Care (PHC), Provincial Health Services Authority (PHSA) and Vancouver Coastal Health (VCH). PHC, PHSA and VCH accept no responsibility for use of this material by any person or organization not associated with PHC, PHSA and VCH. A printed copy of this document may not reflect the current electronic version.