

	RESPIRATORY SERVICES	DATE CREATED: January 2010 DATE REVIEWED/REVISED: October 2015
PROCEDURE	TITLE: <u>Pulmonary Diagnostics:</u> Eucapnic Voluntary Hyperventilation (EVH) (Respiratory Therapy) NUMBER: B-00-12-12126	RELATED DOCUMENTS:

This material has been prepared solely for use at Providence Health Care (PHC), Provincial Health Services Authority (PHSA) and Vancouver Coastal Health (VCH). PHC, PHSA and VCH accept no responsibility for use of this material by any person or organization not associated with PHC, PHSA and VCH. A printed copy of this document may not reflect the current electronic version.

SITE APPLICABILITY:

ST. PAUL'S HOSPITAL

GENERAL INFORMATION:

In May 2001, the IOC Medical Commission convened a workshop to examine asthma, beta-agonists and the Olympic Games. The workshop concluded that any athlete notifying the IOC of the need to use a beta-2 agonist will now require one or more of the following:

- A positive methacholine challenge test
- Evidence of a 12% or greater improvement in FEV₁ following beta-2 agonist
- Bronchial provocation with a Eucapnic Voluntary Hyperventilation (EVH) test
-

After an EVH test, airway hyper-responsiveness (AHR) is confirmed if there is a fall of 10% or more in FEV₁ in the post challenge period. (*Anderson et al 2001, Br J Sports Med 35:344-347*).

EQUIPMENT:

- SensorMedics Vmax Metabolic Cart
- Hans-Rudolf Valve with two silicone connectors
- Silicone mouthpiece
- Microgard Filter
- Noseclip
- Reservoir bag with hose
- H-cylinder of Gas Mixture 5% CO₂, 21% O₂ and balance Nitrogen
- EVH worksheet
- Bronchodilator MDI with spacer

PROCEDURE:

1. Power on the Vmax Metabolic cart, and laptop computer.
2. Open both calibration tanks then calibrate the flow sensor and gas analyzers.
3. From the main menu, choose **New Study** to enter patient information (new patient) or choose **Find Patient** to search for patient information (old patient). Press **F3** to save.
4. From the main menu, choose **Pulmonary Function**, and have patient perform Flow Volume Loops as per ATS guidelines. (Refer to procedure RTD7300). Choose the best FEV₁ and press **F3** to save. The best FEV₁ will be used to calculate ventilation targets and threshold values (EVH worksheet).

5. Fill gas reservoir bag with 5% Carbon Dioxide mixture.
6. On a small table, prepare the Hans-Rudolph valve, mouthpiece and flow sensor as per **Exhibit A**. The Vmax flow sensor will connect to the exhalation side of the Hans-Rudolph valve to allow monitoring of the patient's spontaneous breathing parameters.
7. Prepare the patient to **hyperventilate** at a target range of $(25 - 30) \times FEV_1 = VE$ for **6 minutes**.
8. From the main menu, choose **Exercise/Metabolic Test**. Once in the Metabolic Study screen, choose **F10** to enter the Metabolic Protocol Setup screen.
9. From the Metabolic Protocol Setup screen under the heading **Bike**, choose **External Treadmill**. Under the heading **Text**, choose **EVH Test**. Press **F3** to save changes.
10. Press **Start Test** from the Metabolic Study screen and place the mouthpiece into patient's mouth with nose clips on.
11. Press **F3** again to begin data collection.
12. The patient's breathing should activate the scrolling of the tabular results. In the upper left-hand corner of the screen, click on **Events** and choose **Tech Notes**. In the lower right-hand corner of the screen, choose the **Breath Reject** box to remove the checkmark. This will help decrease the number of breaths rejected by the Vmax during testing. Press **F3** to save changes.
13. Have the patient begin hyperventilation as you toggle into exercise phase at the bottom right-hand corner of the screen. This will begin a counter at the top of the screen. The patient should hyperventilate for 6 minutes. Monitor patient's VE (minute ventilation) to keep it within the calculated target range, giving verbal encouragement as necessary.
14. Adjust flowrate of gas mixture from the H-cylinder to the reservoir bag to keep the bag inflated throughout the 6 minutes of hyperventilation, taking care not to under-fill or over-fill the bag.
15. After 6 minutes of hyperventilation, re-enter **Pulmonary Function** to perform spirometry, and toggle into **Post test**. Perform FEV₁ maneuver at 1, 3, 5, 7, 10, 15, and 20 minutes post hyperventilation.
16. Complete the EVH worksheet and calculate the total percent change achieved.
17. Administer a bronchodilator and re-test until patient returns to their baseline spirometry.



EXHIBIT A

REFERENCES:

1. S D Anderson, G J Argyros, H Magnussen, *et al*. Provocation by eucapnic voluntary hyperpnoea to identify exercise induced bronchconstriction. *British Journal of Sports Medicine* 2001 35: 344-347

REVIEWED BY:

1. Respiratory Therapist, Pulmonary Diagnostics, PHC
2. Pulmonary Diagnostics Coordinator, Respiratory Services, PHC