





ABCD-21-06-90178

Contrast Warmer Temperature Log

Temperature Range of Warmer: 20 degrees C to 30 degrees C.

Alarm to Trigger when below 15 degrees C & above 37 degrees C.

Instructions:

- Technologist to check the temperature of the contrast warmer weekly.
- Do not store food or beverages in the warmer.
- Do not store other substances in the warmer unless stated by the manufacture of the substance.
- Notify supervisor and Biomed if there is an issue with the warmer.

Date (DD-MMM-YYYY)	Current Temperature ©	Initial of MRT	Action/ Follow-up (when temperature is out of range OR warmer is not functioning properly)
			у деграну

Effective date: 29-SEP-2022 Review date: 27-SEP-2022 Document Owner: CT RPL Revision date: 27-SEP-2022

Document Number: ABCD-21-06-90178 Approved by: Regional Director of Quality & Informatics, LMMI