

Temporal Bone or Cholesteatoma Protocol (Contrast)

Field Strength: 1.5T or 3T **Coil:** Multi-channel head coil **Position:** Supine, arms down **Orientation:** Head first

Scan Prescription: Subcallosal plane to prescribe or reformat for axial **Coverage:** Axial - Mastoid Tip to top of Petrous Ridge (*see Appendix A*)
Coronal – Entire Petrous Ridge (*see Appendix A*)

Sequence	FOV (pFOV)	Thickness (Gap)	Resolution	Phase	TR	TE	Comments
Ax T1	180 mm	2 mm (0 mm)	256 x 192	A-P			Fat Sat optional
Cor T1	180 mm	2 mm (0 mm)	256 x 192	S-I			Fat Sat optional
Give Contrast – Single dose per site standards							
Ax T2 Posterior Fossa	220 mm	2 mm (0.5 mm)	256 x 224	A-P			Coverage from foramen magnum to top of Tentorium
Ax 3D CISS / FIESTA-C	160 mm	1.4 mm	320 x 320	A-P			Reformat to coronal
Ax or Cor DWI (non-EPI)	180 mm	3 mm	128 x 128				Axial or coronal plane Non-EPI (HASTE or PROPELLER, etc) 2 b-values = 0-10 and 800-1000 s/mm ²
Ax T1 Fat Sat Post Gad	180 mm	2 mm (0 mm)	256 x 192	A-P			
Cor T1 Fat Sat Post Gad	180 mm	2 mm (0 mm)	256 x 192	S-I			

Notes:

1. This 20-30 minute delay from gad administration to post-gad T1w fat sat images will allow for greater differentiation between types of pathology, if present. Gadolinium does not produce a significant shortening of T2, therefore it is not seen on T2w images post gad of the head.
2. Non-echoplanar DWI is the imaging of choice in detecting cholesteatoma primarily due to the lack of air-bone susceptibility artefact and distortion at the temporal bones (*see Appendix B*).
3. Keratin content in cholesteatoma will return high signal intensity on high b-value compared to non-cholesteatomatous soft tissue such as granulation tissue, inflammation, and fluid, which will return lower or no signal on the high b-value (*see Appendix C*).

Temporal Bone or Cholesteatoma Screening Protocol (Non-Contrast)

Field Strength: 1.5T or 3T **Coil:** Multi-channel head coil **Position:** Supine, arms down **Orientation:** Head first

Scan Prescription: Subcallosal plane to prescribe or reformat for axial **Coverage:** Axial - Mastoid Tip to top of Petrous Ridge (*see Appendix A*)
Coronal – Entire Petrous Ridge (*see Appendix A*)


Sequence	FOV (pFOV)	Thickness (Gap)	Resolution	Phase	TR	TE	Comments
Ax T1	180 mm	2 mm (0 mm)	256 x 192	A-P			
Ax T2 Posterior Fossa	220 mm	2 mm (0.5 mm)	256 x 224	A-P			Coverage from foramen magnum to top of Tentorium
Ax 3D CISS / FIESTA-C	160 mm	1.4 mm	320 x 320	A-P			Reformat to coronal
Ax or Cor DWI (non-EPI)	180 mm	3 mm	128 x 128				Axial or Coronal Non-EPI (HASTE or PROPELLER, etc) 2 b-values = 0-10 and 800-1000 s/mm ²

Notes:


1. CT temporal bone should be used as the initial imaging tool. If tissue confirmation is required, it is recommended that gadolinium be administered to differentiate between post-inflammatory changes and granulation tissue.
2. Non-echo planar DWI is the imaging of choice in detecting cholesteatoma primarily due to the lack of air-bone susceptibility artefact and distortion at the temporal bones (*see Appendix B*).
3. Keratin content in cholesteatoma will return high signal intensity on high b-value compared to non-cholesteatomatous soft tissue such as granulation tissue, inflammation, and fluid, which will return lower or no signal on the high b-value (*see Appendix C*).

Appendices

Appendix A: Coverage for Coronal and Axial images.



Appendix B: EPI DWI versus non-EPI DWI.



Source: Lehmann, P., Brochart, C., Page, C. & Saliou, G. (November 2008). 3T MR Imaging of Postoperative Recurrent Middle Ear Cholesteatomas: Value of Periodically Rotated Overlapping Parallel Lines with Enhanced Reconstruction Diffusion-Weighted MR Imaging. *American Journal of Neuroradiology*. 30(2):423-7

Appendix C: MRI findings on four sequences differentiating cholesteatoma from inflammatory/granulomatous lesions.

Tissue Type	T1W	T2W	Gad Enhancement	Restricted Diffusion
Cholesteatoma	Hypointense	Hyperintense	No	Yes
Granulation Tissue	Hypointense	Hyperintense	Yes	No
Scar	Hypointense/Intermediate	Hypointense/Intermediate	No/Delayed	No
Cholesterol Granuloma	Hyperintense	Hyperintense/Intermediate	No	Variable

Source: Juliano AF, Ginat DT & Mooris GM. (October 2013). Imaging Review of the Temporal Bone. *Radiology*. 269(1):17-33

MRI PROTOCOL

ABCD-21-13-90179

Effective Date:	06-OCT-2022			
Posted Date:	06-OCT-2022			
Last Revised:	05-OCT-2022			
Last Reviewed:	05-OCT-2022			
Approved By: (committee or position)	MIEC	MRI PPC	MRI MPL	
	07-SEP-2022	28-SEP-2022	28-SEP-2022	
Owners: (committee or position)	MRI Regional Practice Lead, LMMI			
Revision History:	Version	Date	Description/ Key Changes	Revised By (Name and Position)
	1.0	06-OCT-2022	Initial Release	Samuel Yim, MRI RPL

This material has been prepared solely for use at Fraser Health (FH), Providence Health Care (PHC), Provincial Health Services Authority (PHSA) and Vancouver Coastal Health (VCH). FH, PHC, PHSA and VCH accept no responsibility for use of this material by any person or organization not associated with FH, PHC, PHSA and VCH. A printed copy of this document may not reflect the current electronic version.