


Quick Release Wrist and Ankle Foam Limb Holders – Disposable RESTRAINT

	<p>Indications:</p> <ul style="list-style-type: none"> Patients assessed to be at risk of disrupting life-saving treatment or where picking, pulling, scratching or peeling exacerbates a skin condition, causes self-injury, or compromises wound site integrity.
<p>Supply Details</p> <p>2 disposable cuffs per package</p> <p>ITEM ID: 00082048</p> <p>VENDOR ITEM ID: Bowers/Posey 2532</p> <p>Manufacturers web page (Video and instructions)</p>	<p>Considerations:</p> <ul style="list-style-type: none"> Remove watches or jewellery before applying limb restraints Side-Lying is the preferred positioning when a patient is in limb restraints. Supine or prone position should be used with clinical judgement Attachment to the bed frame must be made using a quick release buckle as recommended by the manufacturer DO NOT use limb restraints on patients with dislocations, fractures, open wounds or an IV site on the affected limb Ensure the restraint is not applied too tightly or too loosely. You should be able to insert one finger under the ankle or wrist restraint Attach restraints to bed frame only