

# TELEMETRY CENTRAL MONITOR SYSTEM CLEANING LOG

Year: \_\_\_\_\_ Month: \_\_\_\_\_

- Protocol:**
- Clean & disinfect each component a minimum of once daily and initial when complete.
  - Use one wipe for cleaning then use a new wipe for disinfecting.
- Product:**
- Screen can **ONLY** be cleaned & disinfected with Kimtech® wipes or PROSAT® wipes.
  - CaviWipes™ or Accel® wipes can be used on the mouse and keyboard.
- Location:**
- Keep this cleaning log in a plastic sleeve near the telemetry monitor system.

Component:	Telemetry Screen	Keyboard	Mouse
Cleaning frequency:	Minimum once daily	Minimum once daily	Minimum once daily
Date			
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
11			
12			
13			
14			
15			
16			
17			
18			
19			
20			
21			
22			
23			
24			
25			
26			
27			
28			
29			
30			
31			